


# Diagnosis Diabetes in Estonia

## What will provide Hb A1c capillary technology ?

Sébastien SAMMUT, PhD  
Sebia, France

**29/08/2019**




# Diabetes a global health problem

**12%** of global health expenditure is  
spent on diabetes  
(\$673 billion)



# What is diabetes?



- **Diabetes is a chronic condition that occurs when the body cannot produce enough insulin or cannot use insulin**


*Adapted from IDF DIABETES ATLAS  
Eighth edition 2017*

Some figures...


## *Diabetes vs other major diseases*




Diabetes




Hemoglobinopathies



Malaria



HIV/AIDS




Tuberculosis

Major diseases prevalence in the world (WHO)

# Some figures...

## *People with diabetes per country*



Adapted from International Diabetes Federation Atlas (2015)

# Some figures

## *Diabetes in Estonia*

**EESTI DIABEEDILIIT**  
Estonian Diabetes Association




Eestis on diabeet diagnoositud umbes 70 000 inimesel, maailmas põeb diabeeti ca 415 miljonit inimest.

- Diabeet (vana nimega suhkruhaigus) on krooniline ainevahetushaigus, mis vajab igapäevast ja pidevat eneseravi.
- Diabeedi puhul ei tooda kõhunääre insuliini piisavalt või üldse mitte või insuliini mõju organismis on puudulik.
- Eestis on ligikaudu 7000 1 tüübi diabeeti põdevat inimest.
- II tüübi diabeeti põeb teadaolevalt ligikaudu 60-65 000 inimest. Lisaks sellele arvatakse, et umbes sama palju inimesi põeb II astme diabeeti enese teadmata.

- About 70 000 people are diagnosed with diabetes in Estonia
- In addition, it is believed that about the same number of people suffer from stage II diabetes **unknowingly**

# Diabetes: a global health problem

## Adults who died from diabetes, HIV/AIDS, tuberculosis, and malaria



**Figure 1.1**  
The major diabetes complications

### Eye disease

Many people with diabetes develop some form of eye disease (retinopathy), which can damage vision or provoke blindness. Persistently high levels of blood glucose are the main cause of retinopathy. The network of blood vessels that supply the retina can become damaged in retinopathy, leading to permanent loss of vision. Retinopathy however, can become quite advanced before it affects vision, and it is therefore essential that people with diabetes have regular eye screenings. If detected early, treatment can be given to prevent blindness. Keeping good control of blood glucose greatly reduces the risk of retinopathy.

### Cardiovascular disease


Cardiovascular disease is the most common cause of death and disability among people with diabetes. The cardiovascular diseases that accompany diabetes include angina, myocardial infarction (heart attack), stroke, peripheral artery disease and congestive heart failure. High blood pressure, high cholesterol, high blood glucose and other risk factors contribute to the increased risk of cardiovascular complications.

### Pregnancy complications

Women with any type of diabetes are at risk of a number of complications during pregnancy, as high glucose levels can affect the development of the foetus. Women with diabetes therefore require careful monitoring before and during pregnancy to minimise the risk of these complications. High blood glucose during pregnancy can lead to changes in the foetus that cause it to gain excess size and weight. This in turn can lead to problems during delivery, injuries to the child and mother, and low blood glucose (hypoglycaemia) in the child after birth. Children who are exposed to high blood glucose in the womb are at higher risk of developing type 2 diabetes later in life<sup>8</sup>.

### Diabetic foot

As well as nerve damage, people with diabetes can experience problems with poor circulation to the feet, as a result of damage to blood vessels. These problems increase the risk of ulceration, infection and amputation. People with diabetes face a risk of amputation that may be more than 25 times greater than that in people without diabetes<sup>9</sup>. With good management however, a large proportion of amputations can be avoided. Even when a person undergoes amputation, the remaining leg – and the person's life – can be saved by good follow-up care from a multidisciplinary foot team<sup>10</sup>. In view of these risks, it is important that people with diabetes examine their feet regularly.



### Oral health

Diabetes can pose a threat to oral health. There is an increased risk of inflammation of the tissue surrounding the tooth (periodontitis) in people with poor glucose control. Periodontitis is a major cause of tooth loss and is associated with an increased risk of cardiovascular disease. Management of periodontitis is very important in people with diabetes because optimal oral hygiene can prevent tooth loss, facilitate a healthy diet and improve glucose control.

### Kidney disease

Kidney disease (nephropathy) is far more common in people with diabetes than in people without diabetes; diabetes is one of the leading causes of chronic kidney disease. The disease is caused by damage to small blood vessels, which can cause the kidneys to be less efficient, or to fail altogether. Maintaining near-normal levels of blood glucose and blood pressure greatly reduces the risk of nephropathy.

### Prevention of complications

Common to all the major complications of diabetes is the fact that they are not inevitable – they can be prevented by good control of blood glucose levels, as well as good control of blood pressure and cholesterol levels. This requires a high level of education of the person with diabetes in managing their condition, as well as access to insulin, oral medications and monitoring equipment. People with diabetes should be supported by a well-educated health work force as well as health systems that provide regular blood tests and eye and foot examinations. The International Diabetes Federation (IDF) works in many places around the globe to provide treatments and services to improve the outcomes for people with diabetes.

### Nerve damage

Nerve damage (neuropathy) also results from prolonged high blood glucose levels. It can affect any nerve in the body. The most common type is peripheral neuropathy, which mainly affects the sensory nerves in the feet. This can lead to pain, tingling, and loss of sensation. This is particularly significant because it can allow injuries to go unnoticed, leading to ulceration, serious infections and in some cases amputations. Neuropathy can also lead to erectile dysfunction, as well as problems with digestion, urination and a number of other functions.





# Diabetes in Estonia

High probability of premature death (30 - 70 yo)  
from NCDs in Estonia among European  
countries

EUROPEAN JOURNAL OF GENERAL PRACTICE, 2018  
VOL. 24, NO. 1, 112-117  
<https://doi.org/10.1080/13814788.2018.1429594>



Taylor & Francis  
Taylor & Francis Group

BACKGROUND PAPER


OPEN ACCESS



Challenges for clinical practice and research in family medicine in reducing  
the risk of chronic diseases. Notes on the EGPRN Spring Conference 2017  
in Riga


Vija Silina<sup>a</sup> and Ruth Kalda<sup>b</sup>

Family Medicine and Public Health, University of




**Figure 1.** Probability (%) of premature death (dying between age 30 and 70 years) from four main NCDs (cardiovascular diseases, cancer, diabetes, chronic lung disease) in EU countries and in the WHO European region, data from 2015 [5].





# How to diagnose diabetes?



# The guidelines

## Table 2.2—Criteria for the diagnosis of diabetes

FPG  $\geq 126$  mg/dL (7.0 mmol/L). Fasting is defined as no caloric intake for at least 8 h.\*

OR

2-h PG  $\geq 200$  mg/dL (11.1 mmol/L) during OGTT. The test should be performed as described by the WHO, using a glucose load containing the equivalent of 75-g anhydrous glucose dissolved in water.\*

OR

A1C  $\geq 6.5\%$  (48 mmol/mol). The test should be performed in a laboratory using a method that is NGSP certified and standardized to the DCCT assay.\*

OR

In a patient with classic symptoms of hyperglycemia or hyperglycemic crisis, a random plasma glucose  $\geq 200$  mg/dL (11.1 mmol/L).

\*In the absence of unequivocal hyperglycemia, results should be confirmed by repeat testing.

# Diabetes Management

European Journal of General Practice, 2010; 16: 85–91

informa  
healthcare

## ORIGINAL ARTICLE

Meeting targets in type 2 diabetes care contributing to good glycaemic control. A cross-sectional study from a primary care setting in Estonia

ANNELI RÄTSEP<sup>1</sup>, RUTH KALDA<sup>1</sup> & MARGUS LEMBER<sup>2</sup>

<sup>1</sup>Department of Family Medicine, University of Tartu, Estonia, and <sup>2</sup>Department of Internal Medicine, University of Tartu, Estonia

In Estonia 39% of the patients reached the recommended target of HbA1c <6.5% only

## Targets of type 2 diabetes care and associated factors 87




Figure 1. Distribution of type 2 diabetes patients meeting the clinical targets recommended in the guidelines.

# The guidelines


## A1C GOALS

For glycemic goals in children, please refer to Section 12 “Children and Adolescents.”  
For glycemic goals in pregnant women, please refer to Section 13 “Management of Diabetes in Pregnancy.”


### Recommendations

- A reasonable A1C goal for many nonpregnant adults is <7% (53 mmol/mol). **A**
- Providers might reasonably suggest more stringent A1C goals (such as <6.5% [48 mmol/mol]) for selected individual patients if this can be achieved without significant hypoglycemia or other adverse effects of treatment (i.e., polypharmacy). Appropriate patients might include those with short duration of diabetes, type 2 diabetes treated with lifestyle or metformin only, long life expectancy, or no significant cardiovascular disease. **C**
- Less stringent A1C goals (such as <8% [64 mmol/mol]) may be appropriate for patients with a history of severe hypoglycemia, limited life expectancy, advanced microvascular or macrovascular complications, extensive comorbid conditions, or long-standing diabetes in whom the goal is difficult to achieve despite diabetes self-management education, appropriate glucose monitoring, and effective doses of multiple glucose-lowering agents including insulin. **B**

Individualized  
A1c target



## Approach to the Management of Hyperglycemia



**Figure 6.1**—Depicted are patient and disease factors used to determine optimal A1C targets. Characteristics and predicaments toward the left justify more stringent efforts to lower A1C; those toward the right suggest less stringent efforts. Adapted with permission from Inzucchi et al. (72).

# Diabetes Management



ORIGINAL ARTICLE

Severe Hypoglycemia Attributable to Intensive Glucose-Lowering Therapy Among US Adults With Diabetes: Population-Based Modeling Study, 2011-2014

Grace K. Mahoney, MS; Henry J. Henk, PhD; and Rozalina G. McCoy, MD, MS

Individualized  
A1c target will avoid  
hospitalizations because  
of hypoglycemia (9578  
hospitalizations attributed  
to intensive treatment)

## Abstract

**Objective:** To estimate the contemporary prevalence of intensive glucose-lowering therapy among US adults with diabetes and model the number of hypoglycemia-related emergency department (ED) visits and hospitalizations that are attributable to such intensive treatment.

**Patients and Methods:** US adults with diabetes and glycated hemoglobin (HbA<sub>1c</sub>) levels less than 7.0% who were included in the National Health and Nutrition Examination Survey (NHANES) between 2011 and 2014. Participants were categorized as clinically complex if 75 years or older or with 2 or more activities of daily living limitations, end-stage renal disease, or 3 or more chronic conditions. Intensive treatment was defined as any glucose-lowering medications with HbA<sub>1c</sub> levels of 5.6% or less or 2 or more with HbA<sub>1c</sub> levels of 5.7% to 6.4%. First, we quantified the proportion of clinically complex and intensively treated individuals in the NHANES population. Then, we modeled the attributable hypoglycemia-related ED visits/hospitalizations over a 2-year period based on published data for event risk.

**Results:** Almost half (48.8% [10,719,057 of 21,980,034]) of US adults with diabetes (representing 10.7 million US adults) had HbA<sub>1c</sub> levels less than 7.0%. Among them, 32.3% (3,466,713 of 10,719,057) were clinically complex, and 21.6% (2,309,556 of 10,719,057) were intensively treated, with no difference by clinical complexity. Over a 2-year period, we estimated 31,511 hospitalizations and 30,954 ED visits for hypoglycemia in this population; of these, 4774 (95% CI, 954-9714) hospitalizations and 4804 (95% CI, 862-9851) ED visits were attributable to intensive treatment.


**Conclusion:** Intensive glucose-lowering therapy, particularly among vulnerable clinically complex adults, is strongly discouraged because it may lead to hypoglycemia. However, intensive treatment was equally prevalent among US adults, irrespective of clinical complexity. Over a 2-year period, an estimated 9578 hospitalizations and ED visits for hypoglycemia could be attributed to intensive diabetes treatment, particularly among clinically complex patients. Patients at risk for hypoglycemia may benefit from treatment deintensification to reduce hypoglycemia risk and treatment burden.



# Diabetes in Estonia

## Avoidable hospital admissions for diabetes can be improved

**FIG. 7.10** Avoidable hospital admissions for diabetes, asthma and chronic obstructive pulmonary disease (COPD) in selected European countries and Estonia, 2015



Source: OECD Health Statistics, 2017.


Note: Rates are not adjusted by health care needs and health risk factors.

# Diabetes in Estonia


Goal of a quality bonus system (2006) :

- increase the quality and effectiveness of preventive services
- improve monitoring of chronic diseases

**FIG. 3.8** Goal achievements for quality bonus system indicators (between 2006–2014)








# The guidelines

- **According to the ADA:**


- Hb A1c testing must be done on a certified method (NGSP, DCCT, IFCC)
- POC are not recommended for diagnosis
- Interpretation of the Hb A1c value must be done according to factors impacting glycemia, such as hemoglobinopathies!

# What is Hb A1c?



# What is Hb A1c?

What it is:



Hb

Hemoglobin...

A

... more precisely the A fraction ...


1

... which is glycated...

C

... with glucose on the N-terminal valine(s) of the Beta chain(s)

What it is not:



Not bilirubin, Not albumin, ...

Not all the Hb fractions, not Hb S or Hb C, ...

Not glycosylated (labile), not carbamylated, ...

Not with fructose, not on the Alpha chain, not elsewhere on the Beta chain, ...

# What is Hb A1c?

- **Definition (IFCC) :**

- **Hemoglobin (Hb A)** that is **irreversibly glycated** at one or both **N-terminal valines** of the **beta chains**


- **Interests:**

- Hb A1c is a measure of the hemoglobin glycation
- Because the RBC have a lifespan of 120 days, Hb A1c reflects mean glycemia of the last 3 months
- Hb A1c gives info on the long-term glycemic control and allows to monitor the risk of developing complications


# What is Hb A1c?




Hb capillary electrophoresis




**Hb A > 96,5%**




**Hb F < 1%**



**Hb A2 < 3,5%**



# What is Hb A1c?



# How to measure Hb A1c?


- Strict IFCC formula:

$$\% \text{ Hb A1c} = \frac{\text{Hb A1c}}{(\text{Hb A1c} + \text{Hb A0})}$$



# How to measure Hb A1c?

- Hb A1c on capillary electrophoresis:




# How to measure Hb A1c?

- Most common laboratory methods:



Immunoassay

$$[\text{Hb A1c}] + [\text{Hb A2c}] / [\text{Total Hb}]$$



Boronate Affinity

$$\text{Total glycated Hb} / \text{Total Hb}$$



Enzymatic

$$[\text{Hb A1c}] + [\text{Hb A2c}] / [\text{Total Hb}]$$

**“Blind” methods:**  
only a Hb A1c value



HPLC

$$\text{Hb A1c Area} / \sum \text{Hb A Area}$$



Capillary Electrophoresis

$$\text{Hb A1c Area} / (\text{Hb A1c Area} + \text{Hb A0 Area})$$

**Separative methods:**  
Hb A1c value + Hb profile




When choosing/using Hb A1c  
method,  
What should we pay attention to?






# What should we pay attention to?





# What should we pay attention to?





# What should we pay attention to?


## *Analytical performances*

- **Accuracy = Only Imprecision (CV)?**

# What should we pay attention to?

## *Analytical performances*

- Accuracy = Imprecision (CV) and Trueness (*bias*)






# What should we pay attention to?

## *Analytical performances*

- Accuracy = Imprecision (CV) and Trueness (*bias*)




All the important “cut-off” are within 1.2 points.  
CV and bias must be “perfect” in this range of measurements!

# What should we pay attention to?

## *Analytical performances*

- **Accuracy = Imprecision (CV) and Trueness (*bias*)**




High bias: low trueness

# What should we pay attention to?

## *Analytical performances*


- **Accuracy = Imprecision (CV) and Trueness (*bias*)**



# What should we pay attention to?

## *Analytical performances*


- **Accuracy = Imprecision (CV) and Trueness (*bias*)**




# What should we pay attention to?

## *Analytical performances*

- **Accurate Hb A1c method = better patient management**



*Risk of retinopathy progression  
(DCCT, diabetes 1995)*




*Risks for development of various complications  
(Skyler JS, Endocrinol Metab Clin North Am, 1996)*

For example, increase by 1 point the Hb A1c (from 7% to 8%) means multiply by 2 the risk to develop complications after 7 years

# What should we pay attention to?

## *How to assess analytical performances?*



- **IFCC certificates given to manufacturers**


- **Example with capillary electrophoresis instrument:**

- Bias = 0.0 mmol/mol and CV = 0.4 %
- Ranks “Gold” (Biological Variation model)
- $\sigma = 25$  (Sigma-metrics model)

*Goals: CV < 2% NGSP or 3% IFCC, Bronze « Rank » (BV model),  $\sigma > 2$  (Sigma-metrics)*

# What should we pay attention to?

## *How to assess analytical performances?*



**Fig. 2. Models applied to 26 manufacturer/instrument means in CAP 2014 GH2-A survey.**

Mean within-manufacturer interlaboratory CV on the x axis; mean manufacturer absolute bias on the y axis. The black star represents the overall mean of all laboratories. The dots (laboratory instruments) and squares (POCT instruments) represent specific manufacturers with colors for analytical principles: green, immunoassays; red, ion-exchange HPLC; yellow, affinity HPLC; blue, capillary electrophoresis; gray, dry chemistry. Abbott Architect c System (A), Abbott Architect i System (B), Axis-Shield Afinion (C), Beckman AU systems (D), Beckman UniCel Dx C Synchron (E), Bio-Rad D10 (F), Bio-Rad Variant II (G), Bio-Rad Variant II Turbo (H), Bio-Rad Variant Turbo 2.0 (I), Roche Cobas c311 (J), Roche Cobas c500 series (K), Roche Cobas Integra 400 (L), Roche Cobas Integra 800 (M), Sebia Capillarys 2 Flex Piercing (N), Siemens Advia Chemistry Systems (O), Siemens DCA 2000/2000+ (P), Siemens DCA Vantage (Q), Siemens Dimension ExL (R), Siemens Dimension RxL (S), Siemens Dimension Vista (T), Siemens Dimension Xpand (U), Tosoh G7 Auto HPLC (V), Tosoh G8 Auto HPLC (W), Trinity Biotech HPLC (X), Trinity Biotech Premier Hb9210 (Y), Ortho Clin Diag Vitros 5,1 FS, 4600, 5600 Chem System (Z). For more details, see online Supplemental Table 1. opt, Optimum; des, desirable; min, minimum.

### Analysis of the CAP GH2-A survey

- Chemistry (IA)
- Boronate Affinity
- HPLC
- Capillary electrophoresis



# What should we pay attention to?

## *How to assess analytical performances?*





Fig. 3. Quality target models applied to 4 manufacturer/instruments means in the present study and in the CAP 2014 GH2-A survey. Mean within-manufacturer interlaboratory CV on the x axis; mean manufacturer absolute bias on the y axis. The gray stars represent the overall mean of all laboratories in Italy (star inscribed IT) and United States (US). The dots represent the means of analytical devices of manufacturers in both countries: Tosoh (red), Roche (green), Sebia (blue), and Bio-Rad (yellow). (For interpretation of the references to color in this figure legend, the reader is referred to the web version of this article.)

### Analysis of the Italian Hb A1c EQA scheme

- Roche (IA)
- Bio-Rad (HPLC)
- Tosoh (HPLC)
- Sebia (capillary electro.)

# What should we pay attention to?

## *How to assess analytical performances?*



### Analysis of a new European EQA scheme

- Roche (IA)
- Menarini/Arkray (HPLC)
- Tosoh (HPLC)
- Sebia (capillary electro.)
- Bio-Rad (HPLC)
- Siemens (IA)

# What should we pay attention to?


## *Analytical interferences*

- Each instrument is IFCC and NGSP certified





# What should we pay attention to?



# What should we pay attention to?


## *Analytical interferences*

- Triglycerides, total hemoglobin, aspirin, etc
- List of interferences on the NGSP website:
  - <http://www.ngsp.org/interf.asp> and <http://www.ngsp.org/factors.asp>

Method (listed in alphabetical order by manufacturer)	Interference (Yes/No)					
	<i>Hb C trait</i>	<i>Hb S trait</i>	<i>Hb E trait</i>	<i>Hb D trait</i>	<i>Elevated HbF</i>	<i>Carb Hb</i>
Sebia Capillarys 2 Flex Piercing	No 13,59	No 13,59	No 13,54,59	No 13,59	No ≤15% 48	No 48,49,50
Sebia Capillarys 3 Tera	No 13,59	No 13,59	No 13,54,59	No 13,59	-	-





# What should we pay attention to?



# What should we pay attention to?

## *Biological/clinical interferences*

- **Presence of any pathological condition that could impact RBC lifespan and thus the Hb A1c production:**
  - Hemoglobin variants (e.g. Hb S)
  - Thalassemias
- **Around 7% of the worldwide population is carrier for an hemoglobinopathies (WHO)**






# What should we pay attention to?

## *Biological/clinical interferences*

- In Estonia, according to the genetic background (Finnish, Russian...),  $\beta$ -thalassemia, and Beta Hb Variant are found (HbS, Hb E, HbD, HbC)

### Beta Thalassemia

Change in population figure, 2008–2017



In calculating the population figure, Statistics Estonia uses rules for identifying permanent residents, i.e. the residency index. See how register data are used to determine who is a permanent resident and who has left Estonia. More detailed information on the residency index is available in Quarterly Bulletin of Statistics Estonia 1/2017.

About HbVar | Summaries of mutation categories | Query form | Query history | Help | FAQ | Contact us

### HbVar: A database of Human Hemoglobin Variants and Thalassemias

**Query Results**  
There are 21 matches to your query  
Query description: ethnic background in (Russian)

Name	Variant
HbA>NYU	delta 12(A9) Asn>Lys
Initiation codon ATG->ACG (beta)	beta Initiation codon Met>Thr
Codon 15 (G->A): TGG(Trp)->TGA(stop codon) beta0	beta 15(A12) Trp>Stop
Hb Moscow	beta 24(D6) Gly>Asp
Hb Volga	beta 27(B9) Ala>Asp
IVS-I-1 (G->A): AG' GTTGGT->AGATTGGT beta0	beta nt 143 G>A
IVS-I-2 (T->C): AG' GTTGGT->AGACTGGT beta0	beta nt 144 T>C
IVS-I-6 (T->C): the Portuguese type beta+	beta nt 148 T>C
IVS-I-110 (G->A) beta+: the mutation is 21 nucleotides 5' to the acceptor splice site AG'GC	beta nt 252 G>A
IVS-I-130 (G->C): TTAG' GCTG->TTAC' GCTG beta0	beta nt 272 G>C
Hb Tacoma	beta 30(B12) Arg>Ser
Hb Alaska	beta 67(E11) Val>Met
Hb Newcastle	beta 92(F8) His>Pro
Hb Moshaik	beta 92(F8) His>Arg
IVS-II-1 (G->A): beta0	beta nt 496 G>A
IVS-II-654 (C->T): AAGGCAATA->AAG' GTAATA beta+ (reverse)	beta nt 1149 C>T
Hb Durham-NC	beta 114(G16) Leu>Pro
Hb Genoa	beta 139(H17) Asn>Asp
Codon 124-126 (-CCA): Dominant inclusion body: beta-thal trait	Pro: inserted between codons 125(?)
Codon 8 (-AA): AAG(Lys)->-G beta0	beta 8 (-AA): modified C-terminal
Codon 124 (-A): CCA(Pro)->CC: Dominant inclusion body: beta-thal trait	beta 124 (-A): modified C-terminal

<http://globin.bx.psu.edu/hbvar/main.html> August 22, 2019

[http://globin.bx.psu.edu/cgi-bin/hbvar/query\\_vars3](http://globin.bx.psu.edu/cgi-bin/hbvar/query_vars3)

- In Estonia, there is a change in population with a positive net migration



# What should we pay attention to?

## *Biological/clinical interferences*

Blind methods  
(IA, BA, ENZ)

5.1%  
32 mmol/mol

Separative methods  
(CE, HPLC)

5.1%  
32 mmol/mol

# What should we pay attention to?

## *Biological/clinical interferences*

Blind methods  
(IA, BA, ENZ)

5.1%  
32 mmol/mol

Separative methods  
(CE, HPLC)

5.1%  
32 mmol/mol

Just a number. No added value.

# What should we pay attention to?

## *Biological/clinical interferences*


Blind methods  
(IA, BA, ENZ)

5.1%  
32 mmol/mol

Just a number. No added value.

Separative methods  
(CE, HPLC)





5.1%  
32 mmol/mol



Appropriate comment for the clinician

# What should we pay attention to?

## *Biological/clinical interferences*



# What should we pay attention to?

## *Biological/clinical interferences*

### **Effect of RBC Lifespan on HbA1c**

1. Difficult to accurately measure rbc lifespan
2. rbc lifespan ranges 100 - 140 days (average 120 d)
3. Assume rbc lifespan is  $120 \pm 10$  d
4. If HbA1c is 7.0% with lifespan of 120 days
  - 110 days = 6.4%
  - 130 days = 7.6%

*Adapted from Sacks D, AACCC Webinar 2012*

# What should we pay attention to?

## *Biological/clinical interferences*

Tabela

**Análise descritiva de dados de HbA1c, área de Hb, % de área de Hb e glicemia de jejum associados ao perfil cromatográfico por HPLC dos clientes com genótipos de Hb AA, AS e AC obtidas das amostras de sangue total periférico dos 150 clientes participantes do estudo**

Variáveis	(n)	Média	Mediana	Desvio padrão	EPM
A1c HbAA	78	8,05	7,85	1,92	0,21
Área HbAA	–	1768592,37	1729388	396034	44842
GJ HbAA	–	159	147	65,20	7,383
A1c HbAS	40	7,96	7,30	2,18	0,21
Área HbAS	–	2049242,17	2028470	4829	44842
GJ HbAS	–	155	132	80,1	7,383
% Área S	–	34,51	35,05	2,73	0,43
A1c HbAC	32	7,65	7,15	2,26	0,40
Área HbAC	–	2088948,25	1962029	772525	136564
GJ HbAC	–	153	133	81,46	15,97
% Área C	–	30,36	30,45	3,54	0,62

**Patients with Hb A/S or Hb A/C have a lower Hb**

Hb: hemoglobina; HPLC: high-performance liquid chromatography; EPM: erro padrão da média; GJ: glicemia de jejum.



# What should we pay attention to?


## *Biological/clinical interferences*

**Male, 72y**

Presence of abnormal Hb peak = Hb S



**Sickle Cell Trait**, possibly leading to spuriously low Hb A1c value



*Detection of Hb S  
while Hb A1c is measured*

# What should we pay attention to?

## *Biological/clinical interferences*

Table 1. Hemoglobin (Hb) A<sub>1c</sub> (Hb A<sub>1c</sub>) Values Obtained by 4 Different Methods in the Presence of Homozygous and Compound Heterozygous Hb Variants

Hb Variant	Hb A <sub>1c</sub> %			
	Variant II Turbo 1.0 <sup>a,b</sup>	Synchron System <sup>c</sup>	Primus <sup>d</sup>	Integra <sup>e</sup>
SC	24.7	4.7	4.9	6.3
SC	22.0	4.3	4.7	6.1
SC	24.2	4.3	4.1	5.6
SC	26.9	4.6	4.5	6.0
SC	25.9	4.6	4.3	5.5
SS	14.4	4.6	3.3	4.0
SS	13.8	4.4	3.1	3.9
SS	15.8	5.1	4.0	4.5
SS	17.7	4.4	3.1	4.0
SS	23.2	4.2	3.9	4.7
S-β-thalassemia	13.3	4.6	4.9	4.2
S-β-thalassemia	33.0	5.3	5.3	5.8
S-β-thalassemia	10.8	5.4	6.7	7.3
S-β-thalassemia	8.9	5.1	5.4	6.3
S-β-thalassemia	36.3	4.5	4.7	5.4

Patients with Sickle Cell Disease have no Hb A and thus, no Hb A<sub>1c</sub>. But some HPLC, IA and BA methods give a Hb A<sub>1c</sub> results, without any warning!


# What should we pay attention to?

## *Biological/clinical interferences*

### Foetal Hb

- **Increased level of HbF in case of**

- Hereditary persistence of HbF
- $\beta$ -thalassemia
- Sickle cell disease
- Pregnancy
- Anemia
- Leukemias
- Diabetes (Insulin treatment)





# Analytical interference on HPLC




Peak Name	Calibr.	%	Retention Time (min)	Peak Area
A1a	—	1.3	0.170	24718
A1b	—	2.5	0.245	64268
LA1c	—	1.7	0.440	42265
A1c	9.6*	—	0.529	190731
P3	—	3.8	0.786	96362
P4	—	1.5	0.882	37068
Ao	—	82.0	1.044	2069482

$$\% \text{HbA1c} = [100 \times (\text{HbA1c area} / \sum \text{HbA area})]$$

$$\sum \text{HbA area} = (\text{HbA1a}_1/a_2/a_3 \text{ HbA1b} + \text{HbA1c} + \text{P3} + \text{HbA0})$$

A1a<sub>1</sub>+A1a<sub>2</sub> & A1b are minor glycated HbA in which the HbF co-elutes





# IMMUNOASSAYS



$$\%HbA_{1c} =$$

$$[HbA_{1c}] / [Total\ Hb]$$

Total Hb will be falsely increased  
by elevated HbF

Hb A


①

[Total Hb]

Hb F


②

Hb A<sub>2</sub> [Hb A<sub>1c</sub>]






# CAPILLARY ELECTROPHORESIS



Profile with HbF  $\leq 23\%$ , Hb A1c can be reported with a comment after evaluating the clinical context



# What should we pay attention to?

## *Biological/clinical interferences*

- **Presence of any pathological condition that could impact Hb A1c production (clinical and/or analytical interferences):**

- Iron Deficiency Anemia
- Vitamine B12 deficiency
- Antiretroviral treatment
- Hyperthyroidism
- Etc

*Adapted from Gallagher EJ et al,  
Journal of Diabetes 2009*

Factor influencing A <sub>1c</sub>	Increased A <sub>1c</sub>	Decreased A <sub>1c</sub>	Variable change in A <sub>1c</sub>
1. Erythropoiesis	Iron deficiency, vitamin B <sub>12</sub> deficiency, decreased erythropoiesis	Administration of erythropoietin, iron or vitamin B <sub>12</sub> , reticulocytosis, chronic liver disease	
2. Altered hemoglobin			Fetal hemoglobin, hemoglobinopathies, methemoglobin
3. Glycation	Alcoholism, chronic renal failure, decreased erythrocyte pH	Ingestion of aspirin, vitamin C, vitamin E; certain hemoglobinopathies, increased erythrocyte pH	Genetic determinants
4. Erythrocyte destruction	Increased erythrocyte lifespan: splenectomy	Decreased erythrocyte lifespan: hemoglobinopathies, splenomegaly, rheumatoid arthritis, drugs such as antiretrovirals, ribavirin, and dapsone	
5. Assays	Hyperbilirubinemia, carbamylated hemoglobin, alcoholism, large doses of aspirin, chronic opiate use	Hypertriglyceridemia	Hemoglobinopathies

**Figure 1** Factors influencing A<sub>1c</sub>.



# What should we pay attention to?

## *Biological/clinical interferences*

Pathology	Effect on Hb A1c	Effect on Hb
Iron Deficiency Anemia	Spuriously high Hb A1c values	Lowered Hb A2
Vitamine B12 deficiency	Spuriously high Hb A1c values	Elevated Hb A2
Antiretroviral treatment	Spuriously low Hb A1c values	Elevated Hb A2
Hyperthyroidism	Spuriously high Hb A1c values	Elevated Hb A2




Hb A2 can be seen  
and quantified by CE  
and few HPLCs

# What should we pay attention to?

## *Biological/clinical interferences*

### HPLC using short elution program (HbA1c mode)



# What should we pay attention to?

## *Biological/clinical interferences*

**Female, 57y**

MCV = 77.6 fl (!)


Iron = 6.06  $\mu\text{mol/l}$  (!)

Ferritin = 13.3  $\mu\text{g/l}$  (!)

Transferrine = 3.38 g/l (!)



**IDA**, possibly leading to  
spuriously high Hb A1c  
value



*Detection and quantification of the Hb A2  
while Hb A1c is measured*

# Presence of SEBIA Hb A1c capillary technology in the Nordic countries




**(Dark-Green)** Estonia. **(Light-Green)** The rest of the **European Union (EU)**. **(Dark-grey)** The rest of **Europe**. **(Light-grey)** The surrounding region.

# SEBIA Hb A1c capillary technology experience worldwide

- ➡ Allows autonomous overnight running for big laboratories
- ➡ Faster TAT compared to track solution
- ➡ Cost saving solution compared to total laboratory automation
- ➡ Many possibilities of workflow due to multiple techniques on board and the system flexibility

*St Giovanni hospital in  
Florence, Italy  
Two MC3 configurations for  
HbA1c/Hb and P6/IT*



# Sebia support Scientific Activities



***We value and support good scientific practices***



## Main activities

- Technology evaluations
- Laboratory accreditation
- Clinicians' education programs
- User's assistance for interpretation
- Scientific Advisory Board
- KOL collaboration





# The Sebia support

## Customer support – STC and Sebia Academy



- Fully dedicated to total customer satisfaction offering daily assistance for technical and scientific matters
- Second Level support composed by 6 PhDs, highly qualified to answer to the customers needs

### ***SCIENTIFIC AND TECHNICAL SUPPORT***


The Scientific Support Department based in our local Headquarters and at the Sebia affiliates level is dedicated to ensure customer support – 6 Highly qualified PhD

### ***Sebia Academy***

In the last 3 years, the Scientific & Technical Support Department at the Headquarters in France, delivered training to about 5600 people (customers and distributors), both on site and during local events.



# Unique Sebia Hb A1c Separative Solution: We do it right, better, and faster!



Accurate  
and Fast  
Diagnostic

Optimal  
treatment

Less  
complication



# Conclusions




# Conclusion

- **Hb A1c is an important diabetes marker**
  - Long-term glycemia
- **Hb A1c responds to a strict IFCC definition and calculation**
  - Hb A1c is not glycated hemoglobin
- **Hb A1c production is sensible to pathologies affecting RBC lifespan**

# Conclusion

- **When choosing/using Hb A1c method, we should pay attention to:**
  - The analytical performances (CV and bias)
  - The analytical interferences (all conditions that interfere with the measure)
  - The biological/clinical interferences (all condition that can lead to misinterpretation of the Hb A1c value)
- **Only separative technique (such as capillary electrophoresis) bring medical added values**
  - Allow simultaneous detection of several pathological conditions
  - Permit a correct Hb A1c interpretation

Any questions?





Täna tähelepanu eest!

