

Rethinking Platelet Aggregation

Integrated, standardized, available.




Table of contents

Click to jump directly
into the chapter when in
presentation mode...



What are platelets?

3



How do our solutions overcome these challenges?

27



Why is platelet testing so important?

6



What are the advantages for your laboratory?

39



What's behind light transmission aggregation?

16



What makes platelet testing challenging?


24




What are platelets?

What are platelets?

- Smallest blood cells: size 1–2 μm in diameter
- Shaped like a lense (inactive)
- Produced in bone marrow by megakaryocytes
- Lifespan about 5–10 days
- Normal platelet count ranges from 150,000 to 450,000/ μL
- Degraded in liver and spleen by macrophages
- Functions: primary hemostasis, surface for secondary hemostasis, clot retraction



The role of platelets in primary hemostasis





Why is platelet testing so important?

Thrombocytopathy is the most common cause of bleeding

throm·bo·cy·top·a·thy (throm'bō-sī-top'ă-thē)

General term for any disorder of the coagulating mechanism that results from dysfunction of the blood platelets.

Medical Dictionary for the Health Professions and Nursing © Farlex 2012

Inherited vs. acquired thrombocytopathy

Inherited



- Much less frequent
- Can lead to severe bleeding
- Often not diagnosed
- Most common: von Willebrand disease
- Other rare disorders include Glanzmann thrombasthenia, Bernard-Soulier syndrome, and Delta (δ)-storage pool disease

Acquired




- More frequent
- Most common drugs are aspirin, nonsteroidal anti-inflammatory drugs (NSAIDs), along with antiplatelet drugs
- Diseases including cirrhosis, multiple myeloma, kidney disease, and lupus

7% of acutely hospitalized patients show thrombocytopenia at admission^{1,2}

Possible symptoms¹

- Severe bleeding episodes after surgery or injuries
- Hematoma/epitaxis
- Menorrhagia
- Bleeding in the skin/petechiae




1. Scharf RE. Acquired platelet function disorders: pathogenesis, classification, frequency, diagnosis, clinical management. Hämostaseologie. 2008;28:299-311.

2. Guillaume M, et al. Prevalence of thrombocytopenia and thrombocytosis upon acute hospital admission to internal medicine units. A cross-sectional study in Denmark. European Journal of Inner Medicine. 57:e34-e37.

Relation of platelet count and bleeding risk

| Platelet Count | Risk of Bleeding |
|-------------------------------------|--|
| $\geq 50,000/\mu\text{L}$ | Minimal |
| $20,000\text{--}50,000/\mu\text{L}$ | Minor bleeding after trauma |
| $<20,000/\mu\text{L}$ | Spontaneous bleeding |
| $<5000/\mu\text{L}$ | Severe, possibly life-threatening spontaneous bleeding |

→ Reduced platelet function (e.g., due to uremia or aspirin use) adds risk of bleeding in each platelet count range.



Suspect a bleeding disorder?

As a possible diagnostic strategy, ask these important questions:



Is it pathological bleeding, or is the bleeding within the upper normal range?



Is the bleeding due to a congenital/familial or acquired disorder?



Which part of blood coagulation is affected: primary hemostasis (platelet or blood vessel wall problem) or secondary hemostasis (coagulation problem)?




Does the patient have a systemic disease that could cause or exacerbate the bleeding?




Does the patient take any medication that might exacerbate the bleeding tendency?

Steps in bleeding risk screening



Platelet testing scenario



*Not available for sale in the U.S. Product availability varies by country.

Different methodologies for assessment of platelet function

Bleeding time

- 🔴 IVY method
- 🔴 Duke method

Tests based on platelet aggregation under shear stress

- 🔴 INNOVANCE PFA-200*
- 🔴 Impact [cone and platelet analyzer]
- 🔴 Global Thrombosis Test (GTT)

Tests based on platelet aggregation

- 🔴 Light transmission platelet aggregation (LTA)
- 🔴 Impedance aggregation
- 🔴 Lumiaggregometry platelet works

Platelet function methods combined with viscoelastic test

- 🔴 TEG PLATELET MAPPING systems
- 🔴 ROTEM platelets

- 🔴 Uses platelet-rich plasma
- 🔴 Uses whole blood

*Not available for sale in the U.S. Product availability varies by country.
Note: Flow cytometry and thromboxane metabolites are not listed here.



What's behind light transmission aggregometry (LTA)?

Light transmission aggregation: a look back

Gustav Victor Rudolf Born

Born 29 July 1921 in Göttingen

Died 16 April 2018



Gold standard

in platelet function testing




>3300 articles

on light transmission aggregation (LTA)

in PubMed



Light transmission aggregation measurement principle



Source for pictures: Thomas Scholz [2018].

Which reagent is used for which testing purpose?

Differential aggregation responses in patients with thrombocytopathies

| Syndrome | ADP | Collagen | Arachidonic acid | Ristocetin |
|--------------------------|---------|----------|------------------|------------|
| Glanzmann thrombasthenia | ▼ - ∅ | ▼ - ∅ | ▼ - ∅ | n - ▼ |
| Bernard-Soulier syndrome | n | n | n | ▼ - ∅ |
| May-Hegglin syndrome | n - (▼) | n | ? | n |
| δ-Storage pool disease | n - ▼ | ▼ | n - ▼ | n - ▼ |
| α-Storage pool disease | n - ▼ | ▼ | n | n |
| Aspirin-like defect | ▼ | ▼ | ▼ | n |

n = Normal; ∅ = absent; ▼ = inhibited or reduced; (▼) = slightly inhibited.

Disclaimer: When a pathological value is measured in some of the screening tests, further investigations are necessary because different platelet functions (e.g., primary activation, adhesion and aggregation capability, platelet secretion, and procoagulant activity of platelets) may be affected in thrombocytopathy.

Available guidelines for platelet aggregation

1. CLSI. Platelet function testing by aggregometry; approved guideline. CLSI document H58-A. Wayne, PA: Clinical Laboratory Standards Institute; 2008. **[CLSI Guideline]**
2. Cattaneo M, Cerletti C, Harrison P, Hayward CPM, Kenny D, Nugent D, Nurden P, Rao AK, Schmaier AH, Watson SP, Lussana F, Pugliano MT, Michelson AD. Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH. J Thromb Haemost. 2013;11: 1183-9. **[SSC Guideline]**
3. Gesellschaft für Thrombose- und Hämostaseforschung (GTH e.V.): Diagnose von Thrombozytenfunktionsstörungen – Thrombozytopathien, AWMF-S2k-Leitlinie 086-003, 2018. (Available in German language only.) **[GTH Guideline]**


Recommended concentrations according to different guidelines

| | ADP | ARA | COL | EPI | RIS |
|---|---------|------------|-------------|------|-------------------------------|
| SSC 2013 | 2 µM | 1 mM | 2 µg/mL | 5 µM | 1.2 mg/mL / 0.5–0.7 mg/mL |
| CLSI H58-A | 5 µM | 0.5–1.6 mM | 2 µg/mL | 5 µM | 0.8–1.5 mg/mL / ≤0.6 mg/mL |
| German guideline on thrombocytopathies | 2–3 µM* | 1–1.5 mM | 0.5–2 µg/mL | 5 µM | 1.2–1.5 mg/mL / 0.5–0.6 mg/mL |
| Atellica® COAG 360 System[†] | 5 µM | 1 mM | 2 µg/mL | 5 µM | 1.2 mg/mL / 0.5 mg/mL |

*ADP 2–3 µM: Use double concentration if no reaction with low concentration. Evaluation of “second wave”: A sample should aggregate while using 2–3 µM ADP. After holding time, samples require higher concentrations: increase to 4 µM.


[†] Not available for sale in the U.S. Product availability varies by country.

Result reporting: percent aggregation curve



$$\text{Percentage aggregation (t)} = \frac{\text{Start value PRP} - \text{measurement PRP (t)}}{\text{Start value PRP} - \text{PPP value}}$$

Available evaluations for the characterization of the aggregation curve



Result interpretation according to CLSI/ISTH¹

Slope (A), evtl. 2. Slope (A')

Raw value method: Vmax abs, Vmax second wave

Minimum conc. needed to induce secondary wave ADP, EPI (B)

Raw value method: second wave?

Percent maximum aggregation (C)

Raw value method: maximum aggregation

Percent final aggregation (D), disaggregation (C–D)

Raw value method: end aggregation, disaggregation

Lag phase (E)

Raw value method: lag time

Shape change (F)

Raw value method: shape change

Area under the curve (AUC) (G)


Raw value method: AUC

1. CLSI. Platelet function testing by aggregometry; approved guideline. CLSI document H58-A. Wayne, PA: Clinical Laboratory Standards Institute; 2008.



What makes platelet testing challenging?

Challenges in platelet aggregation testing



*“The ideal test of platelet function for use in clinical practice would **be rapid, easy-to-use, inexpensive, and would be a reliable indicator of the clinical response to the specific antiplatelet therapy or combination of therapies.**”*



**How do Siemens Healthineers
integrated platelet aggregation
testing solutions overcome these
challenges?**

Siemens Healthineers and Sysmex integrate platelet aggregation onto routine and specialty hemostasis analyzers



Atellica COAG 360 System is not available for sale in the U.S. Product availability varies by country.



Comprehensive family of systems for platelet aggregation testing

| Hemostasis Portfolio | | | Clotting | Chromogenic | Immunologic | Platelet Aggregation | PSI™ Checks | Automation Connectivity | High-sensitive Immunoassay |
|---|----------------------------|---|---|---|---|---|-------------|-------------------------|---|
| | | |  |  |  |  | | |  |
| High Volume  | Atellica® COAG 360 System* |  | • | • | • | • | • | • | • |
| | Sysmex® CS-5100 System |  | • | • | • | •† | • | • | |
| | Sysmex CS-2500 System |  | • | • | • | •† | • | | |

* Not available for sale in the U.S. Product availability varies by country.

† HYPHEN BioMed CE-marked application. Not available for sale in the U.S. Product availability may vary from country to country and is subject to varying regulatory requirements.

Atellica COAG 360 System*: no more gaps in platelet aggregation testing

Siemens Healthineers platelet aggregation testing solutions can test for five agonists

- HYPHEN BioMed-labeled reagents*
- CE-marked application† and IVDR-conforming
- Available OUS, as RUO in U.S.
- Reagents distributed by Siemens Healthineers and HYPHEN BioMed




Seven Days

7-day dilution of agonists by system and 7-day onboard stability ‡



Ristocetin*




Collagen*,§



Epinephrine*



ADP *



Arachidonic Acid *

* Not available for sale in the U.S. Product availability varies by countries.

† Siemens Healthineers CE-marked applications for Atellica COAG 360 System. Not available for sale in the U.S.

‡ For epinephrine, collagen, and ristocetin.


§ Collagen source: equine reconstituted in acid buffer.

Intelligent sample and reagent management helps to simplify platelet aggregation testing


Atellica COAG 360 System* allows printing of 2D bar codes for reagents.

- The bar-code label enables the system to automatically identify reagent.
- A unique vial number for each reagent ensures monitoring of onboard stability.
- Bar-code maker facilitates the process of final reagent preparation.
- Bar-code labeling can be used for any other reagent.

Sample



Aggregometry



For further instructions, see operator's manual.

*Not available for sale in the U.S. Product availability varies by country.

Available applications for platelet aggregation testing

Siemens Healthineers CE-marked applications

Six CE-marked applications for five different agonists*

- ADP 5 μ M
- Epinephrine 5 μ M
- Arachidonic acid 1 mM
- Ristocetin 1.2 mg/mL and 0.5 mg/mL
- Collagen 2 μ g/mL

Eight additional applications†

- ADP 10 μ M, 2 μ M, 1 μ M
- Epinephrine 10 μ M, 2 μ M, 1 μ M
- Collagen 10 μ M
- TRAP 1 μ M



140 μ L PRP



System requires up to 44% less sample volume than traditional manual methods

*Applications were validated (performance data and application instructions available).

†Without performance data, e.g., reference range.

Atellica COAG 360 System is not available for sale in the U.S.

Improved sample management for platelet aggregation

- Flexible rack configurations: plasma type (PRP/PPP) identified via defined rack
- PPP sample can be run from routine samples
- Advanced pre-analytical sample quality (PSI) checks for PPP samples
- HIL* indices are assigned to the respective aggregation measurement results
- Continuous loading



Three pre-analytical indices

Intelligently identify HIL-affected samples to reduce pre-analytical variables



*HIL: hemolysis, icterus, lipemia


Atellica COAG 360 System is not available for sale in the U.S. Product availability varies by country.

Automated reaction stirring supports more-reliable results in aggregation testing



One standardized process

Automated process streamlines workflow



Four testing channels available on the Atellica COAG 360 System*


*Not available for sale in the U.S. Product availability varies by country.

Integrated platelet aggregation on Atellica COAG 360 System*

Atellica COAG 360 System*




Three manual steps
Operations can be reduced to only three manual steps




*Not available for sale in the U.S. Product availability varies by country.

Atellica COAG 360 System* supports accelerated throughput for platelet aggregation testing

Profile 1: Aggregation only

| | | |
|--------------|-----------|--|
| PRP samples | ● ● ● ● ● |  1 h:25 min 29 PRP samples/h |
| Agonists | ★ ★ ★ ★ ★ | |
| Applications | ■ ■ ■ ■ ■ | |

Profile 2: Typical special hemostasis profile[†] + aggregation

| | | |
|--------------------------------------|---|--|
| Total number of patient samples: 100 | |  2 h:50 min 151 samples/h |
| PRP samples | <div><div></div><div></div><div></div><div></div><div></div></div> | |
| Agonists PRP 1,3,5 | <div><div></div><div></div><div></div><div></div><div></div></div> | |
| Applications PRP 1,3,5 | <div><div></div><div></div><div></div><div></div><div></div><div></div></div> | |
| Agonists PRP 2,4 | <div><div></div><div></div><div></div><div></div><div></div></div> | |
| Applications PRP 2,4 | <div><div></div><div></div><div></div><div></div><div></div><div></div></div> | |

*Not available for sale in the U.S. Product availability varies by country.

[†]PT, aPTT, Fibrinogen, Thrombin, Thrombin Time, INNOVANCE® Antithrombin, INNOVANCE D-Dimer, Protein C, INNOVANCE Free PS Ag, Actin® FS Factor VIII

Atellica COAG 360 System* supports accelerated throughput for platelet aggregation testing

Profile 3: Typical special hemostasis profile[†] + aggregation

Total number of patient samples: 100

| | |
|--------------------|-----------|
| PRP samples | ● ● ● ● ● |
| Agonists PRP 1 | ★ ★ ★ ★ ★ |
| Applications PRP 1 | ■ ■ ■ ■ ■ |
| Agonists PRP 2 | ★ ★ ★ ★ ★ |
| Applications PRP 2 | ■ ■ ■ ■ ■ |



2 h:24 min
163 samples/h

Profile 4: Typical special hemostasis profile[†] + aggregation as STAT

Total number of patient samples: 100

| | |
|--------------------|-----------|
| PRP samples | ● ● ● ● ● |
| Agonists PRP 1 | ★ ★ ★ ★ ★ |
| Applications PRP 1 | ■ ■ ■ ■ ■ |
| Agonists PRP 2 | ★ ★ ★ ★ ★ |
| Applications PRP 2 | ■ ■ ■ ■ ■ |



2 h:06 min
189 samples/h

*Not available for sale in the U.S. Product availability varies by country.

[†]PT, aPTT, Fibrinogen, Thrombin, Thrombin Time, INNOVANCE Antithrombin, INNOVANCE D-Dimer, Protein C, INNOVANCE Free PS Ag, Actin FS Factor VIII

Atellica COAG 360 System* enables thorough traceability of platelet aggregation results

Traceability of results

- Atellica COAG 360 System provides all nine different evaluation parameters for the characterization of the aggregation curve (see slide 18).
- All nine curve parameters can be transferred to the LIS, giving you in-depth audit and archive capabilities.
- Aggregation curve is displayed with customized color coding and comment field.
- Aggregation curves can be printed or exported as PDF via USB,[†] supporting a paperless working environment.



One source

All hemostasis test results come from one source/system



*Not available for sale in the U.S. Product availability varies by country.

† Transfer of PDF files to LIS in development.



What are the advantages for your laboratory?

Siemens Healthineers makes platelet aggregation more manageable, reliable, and accessible

Time- and resource-
consuming testing

No flexible instrument testing

Complex reagent preparation

Variable reproducibility

Large sample volume
required

Result consolidation

Pre-analytical variables

Few specialized labs



*For epinephrine, collagen, and ristocetin.
Atellica COAG 360 System is not available for sale in the U.S.
Product availability varies by country.

Rethink platelet aggregation!
Thank you.



Transforming care delivery through improved access to care.

Atellica, PSI, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. Sysmex is a trademark of Sysmex Corporation. Other trademarks are the property of their respective owners.
The products/features (mentioned herein) are not commercially available in all countries. Due to regulatory reasons their future availability cannot be guaranteed.