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EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance

**Version 1.0
July 2013**

EUCAST subcommittee for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance:

Christian G. Giske (Sweden, EUCAST Steering Committee and EARS-Net Coordination Group; chairman), Luis Martinez-Martinez (Spain, EUCAST Steering Committee), Rafael Cantón (Spain, chairman of EUCAST), Stefania Stefani (Italy), Robert Skov (Denmark, EUCAST Steering Committee), Youri Glupczynski (Belgium), Patrice Nordmann (France), Mandy Wootton (UK), Vivi Miriagou (Greece), Gunnar Skov Simonsen (Norway, EARS-Net Coordination Group), Helena Zemlickova (Czech republic, EARS-Net Coordination Group), James Cohen-Stuart (The Netherlands) and Marek Gniadkowski (Poland).

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Table 1. Clinical breakpoints and screening cut-off values for carbapenemase-producing Enterobacteriaceae (according to EUCAST methodology).


Carbapenem	MIC (mg/L)		Disk diffusion zone diameter (mm) with 10 µg disks	
	S/I breakpoint	Screening cut-off	S/I breakpoint	Screening cut-off
Meropenem ¹	≤2	>0.12	≥22	<25 ²
Imipenem ³	≤2	>1	≥22	<23
Ertapenem ⁴	≤0.5	>0.12	≥25	<25

¹Best balance of sensitivity and specificity

²In some cases zone diameters for OXA-48-producers are up to 26 mm, so **<27 mm** may be used as a screening cut-off during outbreaks caused by OXA-48-producing Enterobacteriaceae, but with reduction in specificity.

Algoritmi projekt

Figure 1. Algorithm for carbapenemase detection.

Abbreviations: APBA=aminophenyl boronic acid, PBA=phenyl boronic acid, DPA=dipicolinic acid (all of them β -lactamase inhibitors added to disks or tablets containing meropenem in combination disk testing assays)

¹ Combination of KPC and MBL can also produce no synergy. Normally these isolates will have very high resistance levels to carbapenems. They are easiest to detect with molecular methods.

² High-level temocillin resistance (MIC >32 mg/L (12, 18), tentative zone diameter ≤ 10 mm with temocillin 30 μ g disk (17)) is a phenotypic indicator of OXA-48 production.

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